

# POWERSTROLL

Wheelchair  
Accessories

**Power**STROLL

*“Lightweight, portable,  
and is available at  
a great price”*



CE

This excellent new Powerstroll will turn most manual wheelchairs into carer controlled powered wheelchairs. It takes the strain out of pushing a heavy wheelchair and makes easy work of hills, ramps and gradients.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

This dual wheel lightweight Powerstroll can easily be fitted and removed within seconds, and adapts to fit most wheelchair sizes.



EASY TO FIT & REMOVE

## Superb Comfort, Style & Performance...

- > Maximum range of 10 miles (16km)\*
- > Fits most wheelchairs with seat widths from 14" to 20" (including narrow track wheelchairs such as the Lightweight Aluminium).
- > Top speed 3mph.\*
- > Fits 2.2cm-2.5cm (7/8"-1") diameter tubing.
- > Uses existing wheelchair handgrips.
- > Lightweight at 14kg and less than 8kg (17lb) with battery pack removed.
- > Twin 'Drive' wheels with 8" black non-marking PU tyres.
- > Removable fabric battery pack.
- > Uses 12v 20ah battery (airline friendly).
- > Lightweight 100w 12 DC motor.
- > Trigger handset with variable speed dial and direction button.
- > Easy to use lift up strap.
- > Weight capacity 115kg (18st).

\*Subject to conditions of use



TRIGGER HANDSET WITH VARIABLE SPEED DIAL & FORWARD/REVERSE SWITCH

CODE: **PWCPP010**

DESCRIPTION: **Lightweight Dual Wheel Powerstroll**



Total Weight:  
14kg



Weight Without Battery:  
7.5kg



Weight Battery Pack Only:  
6.5kg



Max Range:  
10 Miles (16km)\*



Max Gradient:  
6°



Max Weight Capacity:  
115kg (18st)



ADAPTS TO FIT MOST WHEELCHAIRS

Leaflet Code: **LL172**

If you would like to view this product please contact your local stockist, details below.